




# Modified Curved Double Bent Hohmann Retractor

Designed by Lawrence Dorr, MD. Design modification by Bertrand P Kaper, MD

*Designed to be placed on the anterior wall of the acetabulum, the double-bent angle allows for safe retraction of the reflected head of the rectus femoris, reducing concerns of over-retraction*



**PRODUCT NO:**

D6107-MOD

Overall Length: 8.5" (21,6 cm)

Depth from Handle: 4.5" (11,4)

Blade Width: 27.4 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**