



# Modular Head Holder

Designed by Byron E. Dunaway, MD & Wayne Goldstein, MD


*Designed to hold 22mm to 36mm heads for ease of insertion in minimally invasive THR*

7" #8290-01  
9" #8290-02



Head holding ends are plastic coated to help eliminate any damage to the implant. Available in two lengths. Steam and gas sterilizable.

PRODUCT NO'S:
8290-01 [7"] Overall Length: 7" (17,8 cm)
8290-02 [9"] Overall Length: 9" (22,9 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362