



Mueller-type Femoral Neck Elevator with Blunt Teeth

Designed to elevate the proximal femur, the additional blunt teeth on the end allow for better gripping

PRODUCT NO:

3415-01

Overall Length: 13.5" (34,3 cm)

Handle Length: 6.75" (17,1 cm)

Blade Width at Widest: 32 mm



INNOMED 3415-01 CE
(01)00840277128812 (10)1223SS

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362