



## Namba Hip Slide

Designed by Robert S. Namba, MD

*Gently slides the femoral head into the cup without scratching – essential for ceramic heads*

Helps reduce a femoral head trial and implant into the acetabulum during total hip surgery, and helps avoid injury to the knee from excessive traction during hip reduction. Manufactured of delrin to help eliminate damage to the implant. Can be steam or gas sterilized and is radiolucent.

*Facilitates MIS hip replacement procedures*



PRODUCT NO'S:	
Overall Length: 12" (30,5 cm)	
6890	For 22-40 mm heads
6891	For 40-48 mm heads
6892	For 50-60 mm heads

