

Offset Gouge for Posterior Osteophyte Removal in TKA

Designed by Robert Steensen, MD

An offset gouge with a rounded edge designed to more effectively remove osteophytes from round posterior condyles during total knee arthroplasty



PRODUCT NO:

3731

Overall Length: 11" (27,9 cm)
Blade Width: .815" (2 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362