



Designed to remove osteophytes from the posterior femoral condyles during knee arthroplasty

Lotke Offset Osteotome

Designed by Paul Lotke, MD

PRODUCT NO:

4920

Blade Width: 18,5 mm

Overall Length: 9" (22,9 cm)



Wide Offset Osteotome

Designed by Paul Lotke, MD & Adam Rosen, DO

PRODUCT NO:

4935

Blade Width: 13 mm

Overall Length: 9" (22,9 cm)



Dennis Offset Osteotome

Designed by Douglas Dennis, MD & Paul Lotke, MD

PRODUCT NO:

4935-W

Blade Width: 18,5 mm

Overall Length: 9" (22,9 cm)



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362