

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2002 Innomed, Inc., All Rights Reserved.

**Innomed**  
Instrument Division  
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Offset Punches for Total Hip Stem Removal

Two sizes of offsets allow the punches to be used to tap on a distal portion of an intramedullary nail, after a window has been made in the bone below the tip of the nail.

Product No's: 5125-01 [1/2" offset]  
5125-02 [1-1/4" offset]

*Used to help remove an intramedullary nail or prosthesis stem via a window in the shaft of the bone*

