

# Olson Cannulated Swivel Tool


Designed by Daniel Olson, MD

*Designed to help manipulate the capital fragment into corrected alignment for bunion correction after first metatarsal osteotomy*

*Can also be used in joy-stick fashion to manipulate any bone portion/fragment*



Cannulated for use with pins up to .045" (1,1 mm).



**PRODUCT NO:**

**1616-R**

Overall Length: 2.9" (7,4 cm)

T-Handle Width: 1.5" (3,8 cm)

Cannulated for use with Pins: .045" (1,1 mm)

US Patent Pending: D964,565 S



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

[www.innomed.net](http://www.innomed.net)



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362