



OrthoLucent™ Mini Hohmann Retractors

Designed by Jeffrey Lawton, MD

Radiolucent, lightweight retractors

The carbon fiber PEEK material is strong, lightweight, completely radiolucent, can be steam sterilized, and helps to prevent from marring component surfaces.



PRODUCT NO'S:	
1594-R [8mm Blade]	Overall Length: 6.875" (17,5 cm) Blade Width: 8 mm
1597-R [16mm Blade]	Overall Length: 6.875" (17,5 cm) Blade Width: 16 mm

MADE EXCLUSIVELY FOR INNOMED IN SWITZERLAND

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2018 Innomed, Inc.