

Ortho Self-Retaining Retractors

Used to separate the femur and tibia during knee replacement procedures, where the calibrated design can help to balance ligaments

- ▶ Features a no-teeth design, available with flat or serrated outside blades
- ▶ Can also be used for spine surgery where the calibrated ratchet can be used to help accurately measure the size of an opening - useful in procedures to help assess bone graft needs
- ▶ Also useful in foot & ankle surgery

Standard Handle

PRODUCT NO'S:

Flat Outside Pads

1843 [Medium Flat]
Overall Length: 9.25" (23,5 cm)
Blade Width: 10 mm
Blade Thickness: 1.68 mm

1842 [Small Flat]
Overall Length: 6.5" (16,5 cm)
Blade Width: 7 mm
Blade Thickness: 1.68 mm



MADE EXCLUSIVELY
FOR INNOMED IN
G E R M A N Y

PRODUCT NO'S:

Serrated Outside Pads

1843-01 [Medium Serrated]
Overall Length: 9.25" (23,5 cm)
Blade Width: 10 mm
Blade Thickness: 1.68 mm


1842-01 [Small Serrated]
Overall Length: 6.5" (16,5 cm)
Blade Width: 7 mm
Blade Thickness: 1.68 mm



Calibrated ratchet (in mm) is helpful in ligament balancing.



Small Grip Handle



PRODUCT NO'S:

Serrated Outside Pads

1842-01-SG [Small Serrated, Small Grip]
Overall Length: 6.5" (16,5 cm)
Blade Width: 7 mm
Blade Thickness: 1.68 mm



Designed with the grip closer together for easier gripping and to help reduce hand fatigue



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2022 Innomed, Inc.