

# Ortho Self-Retaining Retractors

*The calibrated ratchet can be used to help accurately measure the size of an opening – useful in procedures to help assess bone graft needs*

- ▶ Features a no-teeth design, available with flat or serrated outside blades
- ▶ Useful in knee surgery to separate the femur and tibia during knee replacement procedures, where the calibrated design can help to balance ligaments
- ▶ Also useful in foot & ankle surgery

## Standard Handle

### PRODUCT NO'S:

#### Flat Outside Pads

1843 [Medium Flat]  
Overall Length: 9.25" (23,5 cm)  
Blade Width: 10 mm  
Blade Thickness: 1.68 mm

1842 [Small Flat]  
Overall Length: 6.5" (16,5 cm)  
Blade Width: 7 mm  
Blade Thickness: 1.68 mm



### PRODUCT NO'S:

#### Serrated Outside Pads

1843-01 [Medium Serrated]  
Overall Length: 9.25" (23,5 cm)  
Blade Width: 10 mm  
Blade Thickness: 1.68 mm


1842-01 [Small Serrated]  
Overall Length: 6.5" (16,5 cm)  
Blade Width: 7 mm  
Blade Thickness: 1.68 mm



*Calibrated ratchet (in mm) is helpful assess bone graft needs*



## Small Grip Handle



*Designed with the grip closer together for easier gripping and to help reduce hand fatigue*

### PRODUCT NO'S:

#### Serrated Outside Pads

1842-01-SG [Small Serrated, Small Grip]  
Overall Length: 6.5" (16,5 cm)  
Blade Width: 7 mm  
Blade Thickness: 1.68 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362