


Ortho Suction Tube

Very effective for suction and minor retracting

Helps eliminate plugging due to bone, cement fragments, blood clots, etc.

Designed by
T. Eickmann, MD

Product No: 5465
Overall Length: 9.25" (8.25" from Knurl)
End Hole Dia.: 1mm
Side Hole Dia.: 1.5mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

