

Patella Grasping Forceps

Designed by S. David Stulberg, MD

Bent handle helps the surgeon to evert the patella during minimally invasive knee surgery

Normally two forceps are used. Sold individually.

PRODUCT NO:
4250
Overall Length: 6.75" (17,1 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY




Patella Cover Plate

Designed by S. David Stulberg, MD


Protects the cut surface of the patella during minimally invasive knee surgery

Sharp spikes help hold the plates in place. Lessens the chance of weakening the patella, as pre-drilling is not necessary.

PRODUCT NO'S:
4230-00 [Set of 4 Sizes]
4230-01 [Small] 35 mm x 31 mm
4230-02 [Medium] 36 mm x 32 mm
4230-03 [Large] 37 mm x 33 mm
4230-04 [Extra Large] 38 mm x 34 mm



Flat topside, with three small spikes underneath



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362