

Proximal Tibia Resection Shark Tooth Clamp


Designed by Shara Diers, PA-C

Designed to help grasp and remove the cut proximal portion of the tibia during total and uni knee arthroplasty

PRODUCT NO:

3651
Overall Length: 7" (17,7 cm)
Platform Width: .79" (2 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362