




Quick Clip Rail Clamp

Designed for quick and easy attachment, the clamp snaps on to any OR side rail, and features a pull-to-release lever for quick adjustments



Front



Back

- Compatible with flat-blade style OR accessories and 1" x 1/4" blade post accessories
- Designed for patient weights up to 800 lbs.

PRODUCT NO:
2596
Dimensions:
3" Wide x 2" High x 2.75" Deep (max.)
7,6 cm Wide x 5,1 cm High x 7 cm Deep (max.)



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.