


# OrthoLucent™ Hohmann Retractors

Carbon fiber PEEK composite material is lightweight,  
completely radiolucent and can be steam sterilized

Also helps to prevent from marring component surfaces.



Modified Narrow Hohmann

Modified Blunt Hohmann

Standard Hohmann

Narrow Bent Hohmann

PRODUCT NO'S:
4535-R [Modified Narrow] Blade Width: 18mm Overall Length: 10"
4550-R [Modified Blunt] Blade Width at Widest: 24.5mm Overall Length: 10.75"
4558-R [Standard] Blade Width: 16mm Overall Length: 9.625" (24.4cm)
7110-R [Narrow Bent] Overall Length: 9.75" (23.8cm) Handle Length: 7" (17.8cm) Blade Width: 19mm Depth from Bend: 4.75" (12.1cm)

MADE EXCLUSIVELY  
FOR INNOMED IN  
SWITZERLAND

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362