



Rake Retractor with Contour Grip

Designed for soft tissue retraction, the contour grip allows for a better grip and less fatigue

Features a non-glare finish.



PRODUCT NO:
4839-04 [3-Prong Blade]
Overall Length: 8.125" (20,6 cm)
Rake Width: 13 mm
Rake Depth: 14 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

3-Prong Blade
#4839-04

Rake Retractors with Ergonomic Handle

Designed for soft tissue retraction, the ergonomic handle allows for a better grip and less fatigue

Non-glare finish featured on the metal retractor parts.



PRODUCT NO:
4839 [3-Prong]
Overall Length: 9.5" (24,1 cm)
Rake Width: 13 mm
Rake Depth: 14 mm

USA MADE

3-Prong Blade
#4839

www.innomed.net



INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.