

PROOF




McPherson Retractor Extender

Designed by Ed McPherson, MD

Designed to extend a standard retractor to help provide additional leverage

Now in two sizes to accommodate most retractors – standard for retractors up to .125" (3,2 mm) thick, and large for retractors up to .16" (4 mm) thick.



PRODUCT NO'S:	
6022 [Standard]	Overall Length: 15.625" (39,7 cm) For retractors up to .125" (3,2 mm) thick
6022-01 [Large]	Overall Length: 15.625" (39,7 cm) For retractors up to .16" (4 mm) thick



ISO 13485:2016

Scan to Launch Our Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362