




Rogozinski Glenoid Retractor Set

Designed by Chaim Rogozinski, MD

*Designed with an ergonomic profile to help reduce retraction fatigue and place the assistant's hand out of surgical view, the undersurface of the **Glenoid Retractor** helps stabilize the humeral head to allow excellent visualization of the glenoid, while the **Glenoid Reaming Retractor** helps expose the glenoid for reaming during total shoulder arthroplasty*



Rogozinski Glenoid Retractor
#4271



Rogozinski Glenoid Reaming Retractor
#4277-01

PRODUCT NO:	
4277-00	[Set of One Each]
4271	[Glenoid Retractor]
	Overall Length: 11.25" (28,6 cm)
	Blade Width: .7" (18 mm)
	Prongs Width: 3,4 mm
	Handle Width: 4.5" (11,4 cm)
4277-01	[Glenoid Reaming Retractor]
	Overall Length: 11.25" (28,6 cm)
	Handle Width: 4.5" (11,4 cm)



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362