



Rogozinski Locking Needle Driver/Scissors

Designed by Chaim Rogozinski, MD

Designed with a quick lock & release handle, can drive a needle and cut a suture without changing instruments

PRODUCT NO'S:

3083 [Standard] Overall Length: 6.5" (16,5 cm)
3084 [Large] Overall Length: 7.75" (19,7 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362