

# Roose Utility Knee Retractor

Designed by Paul Roose, DO

*Used for retraction of the soft tissues laterally or medially and for anterior translation of the tibia during tibial prosthetic insertion*

The curvature and width are designed for retraction of soft tissues and excellent visualization of bone structure.



**PRODUCT NO:**

4532

Overall Length: 9"

Blade Width (above tip): 14mm



**FREE TRIAL ON MOST INSTRUMENTS**

© 2011 Innomed, Inc.

**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

**1.800.548.2362**