

"S" Total Knee Retractors


Designed by R. Barry Sorrells, MD

Helps protect the collateral ligaments and popliteal structures while providing excellent visualization within the knee joint

The design is self-retaining and can be used singularly and in pairs. For cruciate sparing or sacrificing prosthetic designs.

PRODUCT NO'S:
3720-00 [Wide Blade] Overall Length: 6" Blade Width: 20mm
3720-01 [Narrow Blade] Overall Length: 6" Blade Width: 10mm


**PROUDLY
MADE
IN THE
USA**



Wide
Blade

*Excellent for minimally
invasive knee procedures*

Narrow
Blade



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362