

# Screw Extractor with Speed Lock

Designed by Khaled Sarraf, MD & Konstantinos Doudoulakis, MD


*Universal extractor designed to accommodate a large range of screws and screw heads from 3.95 to 9.5 mm*

Can also be used to help with removal of other devices that may require a twisting universal locking gripper.



**PRODUCT NO:**

2021  
Overall Length: 9.25" (23,5 cm)  
Jaw Width: 11 mm  
Jaw Length: 5 cm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**