


Silicone Hand with Positioning Rings

Designed to help with positioning of the hand and fingers for surgery, the silicone rings aid in stabilizing the fingers

The flexible silicone is easily bendable, while maintaining the ability to remain in position once set. Silicone hand and rings are steam sterilizable.



PRODUCT NO'S:

1746-00 [Silicone Hand with Pkg of 6 Rings]

Replacement Parts:

1746-02 [Positioning Rings] Package of 12
Dimensions: 28 mm Outer Diameter, 3 mm Thick

MADE FOR INNOMED IN
GERMANY

www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362