



# Stanton Articulating Small Bone Clamps

Designed by John L. Stanton, MD

*Opposing clamps facilitate manipulation of fracture ends*

The small tube allows use of a towel clamp to compress non-union and shortening osteotomies during fixation, as well as to allow the use of Gelpi retractors to distract malunions during revision surgery.




PRODUCT NO'S:	
1811-00	[Set of Left & Right]
<i>Also available individually:</i>	
1811-L	[Left]
Overall Length: 5.125" (13 cm)	
Curved Plate Radius: 5 mm	
Pin Hole for Pins Up To: 2,4 mm	
1811-R	[Right]
Overall Length: 5.125" (13 cm)	
Curved Plate Radius: 5 mm	
Pin Hole for Pins Up To: 2,4 mm	



Towel clamp not included.

*Standard Towel Clamp to Compress  
Gelpi Retractor to Distract*



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Scan to  
Launch Our  
Website



**INNOMED**

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**