


Stoll Retractor/Wire Guide

Designed by Jordan Stoll, MD

Retractor/guide designed for aiming guide wire when performing femoral nailing (TFN for an intertrochanteric fracture), or tibial nailing using the parapatellar approach

PRODUCT NO:

8012
Overall Length: 12" 30,5 cm)
Blade Width: .71" (1,8 cm)
Prong Length: 6 mm
Guide Hole Diameter: 4,5 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362