


Stoll Retractor/Wire Guide

Designed by Jordan Stoll, MD

Retractor/guide designed for aiming guide wire when performing femoral nailing (TFN for an intertrochanteric fracture), or tibial nailing using the parapatellar approach

PRODUCT NO:

8012
Overall Length: 12" 30,5 cm
Blade Width: .71" (1,8 cm)
Prong Length: 6 mm
Guide Hole Diameter: 4,5 mm



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.