



Stulberg Incision Close Gelpi & Blade Set

Designed by S. David Stulberg, MD

Designed to help expose difficult to visualize areas at the end of incisions



PRODUCT NO'S:	
4269-00	[Set - 1 Gelpi & 1 Blade]
Also available Individually:	
4269-01	[Gelpi]
Overall Length: 7.5" (19,1 cm)	
Maximum Spread Width: 3.625" (9,2 cm)	
4269-02	[Blade]
Overall Length: 5.5" (14 cm)	
Blade Width: 1" (2,54 cm)	
Blade Depth from Bend: 1.5" (3,8 cm)	
Blade Bend-Back Angle: 130°	

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED