

Stulberg Proximal Femoral Elevator

Designed by S. David Stulberg, MD

PRODUCT NO:


3420-09

Overall Length: 14" (35,6 cm)

Handle Length: 10" (25,4 cm)

Blade Width at Widest: 48 mm

Blade Width at Prongs: 24 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.