



# Sumko Surgical Finger Guide

Designed by Michael H. Sumko, MD

*Used to help insert a 3.2 mm guide wire, especially during hip fracture surgery, to help prevent puncturing the surgeons' glove*

The entry point for a trochanteric nail can be located through a smaller incision with this device, with reduced risk of penetrating the surgeon's glove while finding the starting point for the guide wire.




US Patent #503638945

**PRODUCT NO:**

8991

Overall Length: 4" (10,2 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**