


# Superior Retractor

Used for retraction around the acetabulum, can be self retaining with the use of 1/8" (3.2 mm) pins



**PRODUCT NO:**  
6027  
Overall Length: 12.75 (32,4 cm)  
Depth from Bend: 5.5" (14 cm)  
Blade Width: 1" (25,4 mm)  
Holes for Pins Up to: 1/8" (3,2 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362