



Sutter Ligament Release Guide

Designed by Dr. med. Damian Sutter

Designed for sonographically guided A1 annular ligament release, the guide is used with a hook knife, and the release is carried out via a small incision while the operation is carried out quasi percutaneously using sonographic control



PRODUCT NO:
1121
Overall Length: 10.1" (25,7 cm)
Insertion Blade Length: 2.125" (5,4 cm)
Insertion Blade Height: 1,7 mm
Insertion Blade Width: 3,4 mm
Guide Channel Width: 2,4 mm



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.