



# Swanson Elevator

Designed by Richard Ferkel, MD

*Angular design helps to go around bone for retraction and elevation – especially useful in small bone surgery of the hand/wrist and foot/ankle*

<b>PRODUCT NO:</b>
1644
Overall Length: 6.375" (16,2 cm)
Blade Depth: .75" (1,9 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

