

T-Handle Femoral Canal Finders

Designed to sound the femoral canal prior to stem broaching, especially useful to help start the broach path during the direct anterior approach

Rockowitz T-Handle Femoral Canal Finder Rasp

Designed by Neal L. Rockowitz, MD

PRODUCT NO:

4990 [Rockowitz Rasp]
Overall Length: 9" (22,9 cm)




ORIGINAL DR. ROCKOWITZ DESIGN Topside Rasp

Rasp on curve topside and sides,
smooth on underside

Modified T-Handle Femoral Canal Finder Rasp

PRODUCT NO:

4989 [Modified Rasp]
Overall Length: 9" (22,9 cm)



MODIFIED DESIGN Underside Rasp

Rasp on curve underside and sides,
smooth on topside

T-Handle Femoral Canal Finder - Smooth

PRODUCT NO:

4990-03 [Smooth]
Overall Length: 9.385" (24,4 cm)



SMOOTH DESIGN

No rasp – smooth underside,
sides, and topside

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362