


# Taper Head Impactor

Designed by Byron E. Dunaway, MD & Wayne Goldstein, MD

*Designed to impact a modular head during minimally invasive THR*



The impactor has a protective coating to interface against the implant to help prevent damage while seating the implant. Can be used with 22 to 32 mm heads. Steam and gas sterilizable.

**PRODUCT NO:**  
7840  
Overall Length: 12" (30,5 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

© 2023 Innomed, Inc.