

Teurlings Medial Malleolar Clamp w/Wire Guide

Designed by Luc Teurlings, MD

Helps to stabilize the medial malleolar fragment during internal fixation



PRODUCT NO:

1803

Cannula Diameter: .062" (1.6mm)
Overall Length: 5.25"

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED The Innomed logo features the word "INNOMED" in a bold, sans-serif font. To the right of the text is a stylized lowercase "i" that incorporates a green and blue circular design. Further to the right is a small American flag icon.

1.800.548.2362