


# Teurlings Modified Bent Hohmann Retractor

Designed by Luc Teurlings, MD

*Designed to help protect the femur cuts while retracting the MCL*

The twisted blunt end also helps elevate the femur and protect the MCL.



**PRODUCT NO:**

7109

Overall Length: 8.625" (21,9 cm)

Depth from Bend: 4.25" (10,8 cm)

Blade Width (Upper Portion): 17 mm

Blade Width (Lower Portion): 11 mm



ISO 13485:2016  
**FREE TRIAL ON MOST INSTRUMENTS**

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**