



Tissue Protector

Helps protect tissue when a straight reamer is being used

Designed to be used when a straight reamer is being used in a bone canal. Very useful in minimally invasive total hip arthroplasty.



PRODUCT NO'S:	
5480-01	Inside Diameter: 19 mm Overall Length: 6.5" (16,5 cm) Tube Depth: 3.875" (9,8 cm)
5480-02	Inside Diameter: 24 mm Overall Length: 6.5" (16,5 cm) Tube Depth: 3.875" (9,8 cm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.