


Torx/Hex Adapter Set

Designed by Stephen M. Walsh, MD

Designed for conversion of a 3.5 mm screwdriver

Especially helpful when an articulated, universal joint driver is needed (i.e. acetabular screws)



Hex Bit to Torx Driver Adapter

Torx Bit to Hex Driver Adapter

PRODUCT NO'S:
8003-00 [Set – One Each]
Set Includes/Available Separately:
8003-01 [Torx Bit to Hex Driver Adapter] Overall Length: .6" (1,54 cm)
8003-02 [Hex Bit to Torx Driver Adapter] Overall Length: .6" (1,54 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362