


# Universal Traction Assembly

A universal traction assembly with Traction Device for standard operating room tables used to assist with fracture fixation in the acetabulum, pelvis, and femur, and designed to attach to standard operating table side rails

<b>PRODUCT NO'S:</b>	
<b>0006-00</b>	[Universal Traction Assembly]
<b>Set Includes / Available Individually:</b>	
<b>0006-01</b>	[Universal Traction Device]
<b>0005-02</b>	[Vertical Rod]
<i>One with this product number, two included in set.</i>	
<b>0005-03</b>	[Horizontal Rod]
<b>0008</b>	[Disposable Sterile Kit] <i>Kit Includes:</i>
<i>(1) Impervious Stockinette and (1) 11 ft. Traction rope</i>	
<b>Optional Items:</b>	
<b>0008-CASE</b>	[Case of Sterile Kits] <i>Pkg of 10</i>
<b>2595</b>	[Table Clamp]



All components can be ordered separately.



[www.innomed.net](http://www.innomed.net)



**FREE TRIAL ON MOST INSTRUMENTS**

**INNOMED**

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**