



# Vaughan Endzone Retractor


Designed by Roderick Vaughan, MD

*Designed for use when placing the end screws while plating a fracture using a minimally invasive technique*

The "U"-shaped wall design helps allow the maximal exposure along the length, or "endzone", of an incision while maintaining adequate width and retraction along the sides of the exposure.

**PRODUCT NO:**

1766  
 Overall Length: 8.75" (22.2 cm)  
 Deep Depth: 45 mm  
 Deep Internal Width: 14 mm  
 Shallow Depth: 25 mm  
 Shallow Internal Width: 12 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
 www.innomed.net info@innomed.net

912.236.0000 Phone  
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
 Fax +41 41 740 67 71

**1.800.548.2362**

© 2019 Innomed, Inc.