


Self-Centering Verbrugge Bone Clamp

*Self centering Verbrugge forceps with
easy release locking mechanism*



PRODUCT NO:
3639
Overall Length: 8" (20,3 cm)
Maximum Jaw Opening: 2" (5,1 cm)
Upper Jaw Width: 7 mm
Lower Jaw Width: 2,5 mm

MADE FOR INNOMED IN
GERMANY

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362