




Weinraub Joint and Calcaneal Spreader

Designed by Glenn M. Weinraub DPM, FACFAS

Designed to assist in the opening of small joints of the foot and hand for the application of fusion and graft techniques

Provides excellent joint exposure without blocking intra-articular or osteotomy access. Helps prevent slippage or falling out of the joint by placing the arms on either side of the area to be distracted, driving two pins and opening the joint.



Speed lock helps allow precise control and prevent unintended release

PRODUCT NO'S:	
1870 [Standard 1.6 mm]	Overall Length: 7" (17,8 cm) Pin Diameter: Up to .062" (1/16") (1.6 mm)
1872 [Standard 2.8 mm]	Overall Length: 7" (17,8 cm) Pin Diameter: Up to .11" (7/64") (2.8 mm)
1870-SL [Speed Lock 1.6 mm]	Overall Length: 7" (17,8 cm) Pin Diameter: Up to .062" (1/16") (1.6 mm)
1872-SL [Speed Lock 2.8 mm]	Overall Length: 7" (17,8 cm) Pin Diameter: Up to .11" (7/64") (2.8 mm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

®