

# Wells Modified Lambotte PAO Osteotomes

Designed by Joel Wells, MD

*Designed to focus on the posterior column osteotomy and connection to the ischial cut – straight, curved and two offset options helps the posterior column osteotomy to be cut with more control*

Silicone handle designed for better control.



PRODUCT NO'S:	
<b>5276-00</b>	[Set of 1 each with case]
<b>Set Includes / Available Individually:</b>	
<b>5276-01</b>	[Straight PAO Osteotome] Overall Length: 14" (35,6 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
<b>5276-02</b>	[Curved PAO Osteotome] Overall Length: 13.875" (35,3 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
<b>5276-03</b>	[Offset PAO Osteotome] Overall Length: 13.625" (34,6 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
<b>5276-04</b>	[Offset Curved PAO Osteotome] Overall Length: 13.625" (34,6 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
<b>9007</b>	[Case]

MADE EXCLUSIVELY  
FOR INNOVATED BY  
GERMANY

[www.innomed.net](http://www.innomed.net)



# INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.