


Wells Modified Lambotte PAO Osteotomes

Designed by Joel Wells, MD

Designed to focus on the posterior column osteotomy and connection to the ischial cut – straight, curved and two offset options helps the posterior column osteotomy to be cut with more control

Silicone handle designed for better control.



PRODUCT NO'S:	
5276-00	[Set of 1 each with case]
Set Includes / Available Individually:	
5276-01	[Straight PAO Osteotome] Overall Length: 14" (35,6 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 11,1 mm
5276-02	[Curved PAO Osteotome] Overall Length: 13.875" (35,3 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
5276-03	[Offset PAO Osteotome] Overall Length: 13.625" (34,6 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
5276-04	[Offset Curved PAO Osteotome] Overall Length: 13.625" (34,6 cm) Handle Length: 4.5" (11,4 cm) Osteotome Width: 12,7 mm
9007	[Case]

MADE EXCLUSIVELY
FOR INNOVATED BY
GERMANY

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362