


# Wiater Shoulder Drape Crossbar

Designed by J. Michael Wiater, MD, FAOS, FAOA


*Designed for use during shoulder surgery in the beach chair position or during other surgical procedures to support and keep the surgical drapes away from the surgical site, maintain a sterile field, and help to allow the anesthesia provider good access to the airway*



**PRODUCT NO:**  
**2417**  
Overall Length: 60" (152,4 cm)



Lightweight stainless steel bar with end clamps for attaching to two IV poles.



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**