



# Wilson Double Scalpel Handle

Designed by Ralph Wilson, MD

*Designed to help make a predictable incision in the patellar tendon when harvesting ACL graft material*

The blade offset is 10 millimeters. The tendon graft is harvested from the patella and tibial tubercle including the patellar tendon. Uses scalpel blades that fit a #3 handle size.



Scalpel blades not included.

**PRODUCT NO:**  
8207  
Overall Length: 5.75" (14,6 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**