


Witmer Deep APC Anterior Hip Retractor

Design modification by Dan Witmer, MD of an original design by APC, Inc.

Designed with extra length and depth for use around the acetabulum with large patients during total hip replacement surgery



PRODUCT NO:

6420-01

Overall Length: 14" (35,6 cm)
Upward Bent Handle Length: 5.45" (13,8 cm)
Handle-to-2nd-Bend Length: 13" (16,5 cm)
Blade Depth: 7.15" (18,1 cm)
Blade Width: 22 mm



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362