

Wubben Lateral Fat Pad Retractor for TKR


Designed by Robert Wubben, MD

Designed to hold soft tissues when inserting the TKR

PRODUCT NO:

3218

Overall Length: 10" (25,4 cm)
Blade Width: 41 mm



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS



ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.