



Wurapa Small Joint Compressor and Distractor

Designed by Raymond K. Wurapa, MD

Designed to allow one-handed manipulation and deployment once fixation pins are placed

Designed to simplify several small joint procedures:

- ▶ Preparation of small bone non-unions before bone grafting and fixation
- ▶ Preparation of small joints for arthrodesis (e.g. partial wrist fusion)
- ▶ Useful for intercarpal stabilization while performing ligament reconstructions (e.g. scapholunate ligament repair/reconstruction)
- ▶ Distract and better evaluate small joints before determining final management.



Available with two hole sizes on each instrument!

Pins should be cut short above the pin guides to allow full access to the operative site.



Distractor with One Hole Size
#1754



Distractor with Two Hole Sizes
#1752



Compressor with One Hole Size
#1753



Compressor with Two Hole Sizes
#1751

PRODUCT NO'S:

DOUBLE HOLES: .045" (1,1 mm) & .062" (1,6 mm)

SINGLE HOLE: .045" (1,1 mm) Hole

1751 [Compressor]
Compresses From: 28 mm
Overall Length: 4.625" (11,7 cm)

1753 [Compressor]
Compresses From: 28 mm
Overall Length: 4.5" (11,4 cm)

1752* [Distractor]
Distracts to: 46 mm
Overall Length: 4.625" (11,7 cm)

1754 [Distractor]
Distracts to: 46 mm
Overall Length: 4.5" (11,4 cm)

MADE EXCLUSIVELY FOR INNOMED IN GERMANY



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362